


Product Datasheet

MCAK (KIF2C) antibody (orb1322934)

Description	Mouse monoclonal antibody to MCAK (KIF2C)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human KIF2C(NP_006836) produced in HEK293T cell.
Target	KIF2C
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	81.1 kDa
Uniprot ID	Q99661
NCBI	NM_006845
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...