


Product Datasheet

PTP1B (PTPN1) antibody (orb1322943)

Description	Mouse monoclonal antibody to PTP1B (PTPN1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PTPN1(NP_002818) produced in HEK293 cell.
Target	PTPN1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	49.8 kDa
Uniprot ID	P18031
NCBI	NM_002827
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...