

Product Datasheet

PTP1B (PTPN1) antibody (orb1322944)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to PTP1B (PTPN1) |
| Species/Host | Mouse |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Full length human recombinant protein of human PTPN1(NP_002818) produced in HEK293T cell. |
| Target | PTPN1 |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| MW | 49.8 kDa |
| Uniprot ID | P18031 |
| NCBI | NM_002827 |
| Dilution Range | WB 1:2000, IHC 1:150 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...