


Product Datasheet

Beclin 1 (BECN1) antibody (orb1323470)

Description	Mouse monoclonal antibody to Beclin 1 (BECN1)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human BECN1 (NP_003757) produced in HEK293T cell.
Target	BECN1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	51.7 kDa
Uniprot ID	Q14457
NCBI	NM_003766
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...