


## Product Datasheet

### Beclin 1 (BECN1) antibody (orb1323490)

<b>Description</b>	Mouse monoclonal antibody to Beclin 1 (BECN1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, WB
<b>Immunogen</b>	Full length human recombinant protein of human BECN1 (NP_003757) produced in HEK293T cell.
<b>Target</b>	BECN1
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	51.7 kDa
<b>Uniprot ID</b>	<a href="#">Q14457</a>
<b>NCBI</b>	<a href="#">NM_003766</a>
<b>Dilution Range</b>	WB 1:500~2000, IF 1:100, FLOW 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...