


Product Datasheet

Aquaporin 1 (AQP1) antibody (orb1323717)

Description	Mouse monoclonal antibody to Aquaporin 1 (AQP1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human AQP1 (NP_932766) produced in HEK293T cell.
Target	AQP1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	28.3 kDa
Uniprot ID	P29972
NCBI	NM_198098
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...