


Product Datasheet

Adenylosuccinate Lyase (ADSL) antibody (orb1324200)

Description	Mouse monoclonal antibody to Adenylosuccinate Lyase (ADSL)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human ADSL (NP_000017) produced in HEK293T cell.
Target	ADSL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54.7 kDa
Uniprot ID	P30566
NCBI	NM_000026
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...