


Product Datasheet

RDH11 antibody (orb1324893)

Description	Mouse monoclonal antibody to RDH11
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human RDH11(NP_057110) produced in HEK293T cell.
Target	RDH11
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	35.2 kDa
Uniprot ID	Q8TC12
NCBI	NM_016026
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...