


Product Datasheet

Transketolase (TKT) antibody (orb1325045)

Description	Mouse monoclonal antibody to Transketolase (TKT)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TKT (NP_001055) produced in HEK293T cell.
Target	TKT
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	67.9 kDa
Uniprot ID	P29401
NCBI	NM_001064
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...