


## Product Datasheet

### Annexin A10 (ANXA10) antibody (orb1325055)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Annexin A10 (ANXA10)                                                                                       |
| <b>Species/Host</b>        | Mouse                                                                                                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | IF, WB                                                                                                                                  |
| <b>Immunogen</b>           | Full length human recombinant protein of human ANXA10(NP_009124) produced in HEK293T cell.                                              |
| <b>Target</b>              | ANXA10                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG2a                                                                                                                                   |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>MW</b>                  | 37.3 kDa                                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9UJ72</a>                                                                                                                  |
| <b>NCBI</b>                | <a href="#">NM_007193</a>                                                                                                               |
| <b>Dilution Range</b>      | WB 1:1000, IF 1:100                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Western blot analysis of HEK293T cells w...