

## Product Datasheet

### GRASP1 (GRIPAP1) antibody (orb1325061)

**Description** Mouse monoclonal antibody to GRASP1 (GRIPAP1)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** FC, IF, IHC, WB

**Immunogen** Full length human recombinant protein of human GRIPAP1 (NP\_064522) produced in HEK293T cell.

**Target** GRIPAP1

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG1

**Clonality** Monoclonal


**MW** 96.0 kDa

**Uniprot ID** [Q4V328](#)

**NCBI** [NM\\_020137](#)

**Dilution Range** WB 1:1000, IHC 1:50, IF 1:100, Flow 1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...