


Product Datasheet

GRASP1 (GRIPAP1) antibody (orb1325078)

Description	Mouse monoclonal antibody to GRASP1 (GRIPAP1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human GRIPAP1 (NP_064522) produced in HEK293T cell.
Target	GRIPAP1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96.0 kDa
Uniprot ID	Q4V328
NCBI	NM_020137
Dilution Range	WB 1:1000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...