


Product Datasheet

B Raf (BRAF) antibody (orb1325081)

Description	Mouse monoclonal antibody to B Raf (BRAF)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	Full length human recombinant protein of human BRAF (NP_004324) produced in HEK293T cell.
Target	BRAF
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	84.4 kDa
Uniprot ID	P15056
NCBI	NM_004333
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...