


Product Datasheet

Peroxiredoxin 4 (PRDX4) antibody (orb1325235)

Description	Mouse monoclonal antibody to Peroxiredoxin 4 (PRDX4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PRDX4 (NP_006397) produced in HEK293T cell.
Target	PRDX4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	30.4 kDa
Uniprot ID	Q13162
NCBI	NM_006406
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...