


Product Datasheet

Cytochrome P450 Reductase (POR) antibody (orb1325283)

Description	Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human POR (NP_000932) produced in HEK293T cell.
Target	POR
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	76.9 kDa
Uniprot ID	P16435
NCBI	NM_000941
Dilution Range	WB 1:1000, IHC 1:50, IF 1:50~100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...