


Product Datasheet

Cytochrome P450 Reductase (POR) antibody (orb1325302)

Description	Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human POR expression vector
Target	POR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	76.7 kDa
Uniprot ID	P16435
NCBI	NM_000941
Dilution Range	WB 1:1000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...