

## Product Datasheet

# Cytochrome P450 Reductase (POR) antibody (orb1325303)

**Description** Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** FC, IF, IHC, WB

**Immunogen** Full length human recombinant protein of human POR (NP\_000932) produced in HEK293T cell.

**Target** POR

**Note** For research use only

**Isotype** IgG2a

**Clonality** Monoclonal


**MW** 76.9 kDa

**Uniprot ID** [P16435](#)

**NCBI** [NM\\_000941](#)

**Dilution Range** WB 1:500, IHC 1:50, IF 1:50, FLOW 1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...