


## Product Datasheet

# Cytochrome P450 Reductase (POR) antibody (orb1325304)

<b>Description</b>	Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human POR (NP_000932) produced in HEK293T cell.
<b>Target</b>	POR
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	76.9 kDa
<b>Uniprot ID</b>	<a href="#">P16435</a>
<b>NCBI</b>	<a href="#">NM_000941</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:50, IF 1:50~100, FLOW 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...