


Product Datasheet

Isocitrate dehydrogenase (IDH1) antibody (orb1325305)

Description	Mouse monoclonal antibody to Isocitrate dehydrogenase (IDH1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human IDH1 expression vector
Target	IDH1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.7 kDa
Uniprot ID	O75874
NCBI	NM_005896
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...