


## Product Datasheet

### OXSR1 antibody (orb1325515)

<b>Description</b>	Mouse monoclonal antibody to OXSR1
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Full-length protein expressed in 293T cell transfected with human OXSR1 expression vector
<b>Target</b>	OXSR1
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	57.8 kDa
<b>Uniprot ID</b>	<a href="#">O95747</a>
<b>NCBI</b>	<a href="#">NM_005109</a>
<b>Dilution Range</b>	WB 1:1000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...