

Product Datasheet

Purified HIV1 p24 antibody (orb1325671)

Description	Mouse monoclonal antibody to Purified HIV1 p24
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-211 of HIV gag p24
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.