


Product Datasheet

SARS-CoV-2 N Protein antibody (orb1325780)

Description	Mouse monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Full length recombinant SARS-CoV-2 N protein produced in HEK293.
Target	SARS-CoV-2 N Protein
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	47kDa
Dilution Range	WB 1:500-1000, ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of full length rec...