

Product Datasheet

S Protein antibody (orb1325795)

Description	Mouse monoclonal antibody to S Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Spike-RBD protein produced in E.coli.
Target	S Protein
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	34kDa
NCBI	YP_009724390
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.