


Product Datasheet

BCMA (TNFRSF17) antibody (orb1325997)

Description	Mouse monoclonal antibody to BCMA (TNFRSF17)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant full length protein of human TNFRSF17/BCMA (NP_001183) produced in HEK293 cell.
Target	TNFRSF17
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.1 kDa
Uniprot ID	Q02223
NCBI	NM_001192
Dilution Range	FC 1:100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T ...