



Product Datasheet

BDKRB2 antibody (orb13265)


Description	BDKRB2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Bdkrb2 (51-150/391 aa)
Target	BDKRB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P30411
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of BDKRB2 antibody staining in cells.



Paraffin-embedded rat brain tissue fixed and stained with BDKRB2 antibody.



Paraffin-embedded rat brain tissue fixed and stained with BDKRB2 antibody.