


Product Datasheet

CD272 (BTLA) antibody (orb1326598)

Description	Mouse monoclonal antibody to CD272 (BTLA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant protein fragment corresponding to amino acids 31-157+Mouse Fc of human BTLA (NP_861445) produced in HEK293T cell.
Target	BTLA
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q7Z6A9
NCBI	NM_181780
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T ...