


Product Datasheet

LDL Receptor (LDLR) antibody (orb1326652)

Description	Mouse monoclonal antibody to LDL Receptor (LDLR)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human LDLR (NP_000518) produced in HEK293T cell.
Target	LDLR
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	95.38 kDa
Uniprot ID	P01130
NCBI	NM_000527
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T

...