



Product Datasheet

B7-DC antibody (orb13269)

Description	B7-DC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human PDCD1LG2 (41-150/273 aa)
Target	PDCD1LG2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	Q9BQ51
Dilution Range	ICC=1:100, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



IHC-P of human cervical cancer tissue (B...



Immunohistochemical analysis of paraffin...