


Product Datasheet

CD229 (LY9) antibody (orb1326968)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody to CD229 (LY9) |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 48-454 of human LY9 (NP_002339) produced in SF9 cell. |
| Target | LY9 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 72 kDa |
| Uniprot ID | Q9HBG7 |
| NCBI | NM_002348 |
| Dilution Range | WB 1:2000, IHC 1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...