

Product Datasheet

ESE1 (ELF3) antibody (orb1328146)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody to ESE1 (ELF3) |
| Species/Host | Mouse |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 175-286 of human ELF3(NP_001107781) produced in E.coli. |
| Target | ELF3 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| MW | 41.3 kDa |
| Uniprot ID | P78545 |
| NCBI | NM_001114309 |
| Dilution Range | WB 1:2000 |
| Expiration Date | 12 months from date of receipt. |

