


Product Datasheet

TAK1 (MAP3K7) antibody (orb1328350)

Description	Mouse monoclonal antibody to TAK1 (MAP3K7)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 254-367 of human MAP3K7(NP_003179) produced in E.coli.
Target	MAP3K7
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	64 kDa
Uniprot ID	O43318
NCBI	NM_003188
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...