


## Product Datasheet

### XPB (ERCC3) antibody (orb1328524)

<b>Description</b>	Mouse monoclonal antibody to XPB (ERCC3)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 191-436 of human ERCC3 (NP_000113) produced in E.coli.
<b>Target</b>	ERCC3
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	89.1 kDa
<b>Uniprot ID</b>	<a href="#">P19447</a>
<b>NCBI</b>	<a href="#">NM_000122</a>
<b>Dilution Range</b>	WB 1:500, IHC 1:1000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of extracts (35ug)...