


Product Datasheet

C20orf7 (NDUFAF5) antibody (orb1328657)

Description	Mouse monoclonal antibody to C20orf7 (NDUFAF5)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human NDUFAF5 (NP_077025) produced in HEK293T cell.
Target	NDUFAF5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	38.7 kDa
Uniprot ID	Q5TEU4
NCBI	NM_024120
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...