

Product Datasheet

12 Lipxygenase (ALOX12) antibody (orb1329535)

Description

Mouse monoclonal antibody to 12

170 —
130 —


Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

IHC

Immunohistochemical staining of paraffin...

Immunogen

Full length human recombinant protein of human ALOX12 (NP_000688) produced in HEK293T cell.

Target

ALOX12

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype

IgG1

Clonality

Monoclonal

MW

75.5 kDa

Uniprot ID

P18054

NCBI

NM_000697

Dilution Range

IHC 1:150

Expiration Date

12 months from date of receipt.