
Product Datasheet

Glucocorticoid Receptor (NR3C1) antibody (orb1330245)

Description

Mouse monoclonal antibody to Glucocorticoid



Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IHC, WB

Western blot analysis of HEK293T cells w...

Immunogen

Human recombinant protein fragment corresponding to amino acids 1-265 of human NR3C1(NP_000167) produced in E.coli.

Target

NR3C1

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype

IgG1

Clonality

Monoclonal

MW

85.5 kDa

Uniprot ID

[P04150](#)

NCBI

[NM_000176](#)

Dilution Range

WB 1:2000, IHC 1:150

Expiration Date

12 months from date of receipt.