


Product Datasheet

Galectin 3 (LGALS3) antibody (orb1330938)

Description	Mouse monoclonal antibody to Galectin 3 (LGALS3)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human LGALS3 (NP_002297) produced in E.coli.
Target	LGALS3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	P17931
NCBI	NM_002306
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...