


Product Datasheet

Dicer (DICER1) antibody (orb1332340)

Description	Mouse monoclonal antibody to Dicer (DICER1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1666-1922 of human DICER1 (NP_085124) produced in E.coli.
Target	DICER1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	218.5 kDa
Uniprot ID	Q9UPY3
NCBI	NM_030621
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HT29 cell lysat...