


Product Datasheet

TRIF (TICAM1) antibody (orb1332762)

Description	Mouse monoclonal antibody to TRIF (TICAM1)
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 120-420 of human TICAM1 (NP_891549) produced in E.coli.
Target	TICAM1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	76.2 kDa
Uniprot ID	Q8IUC6
NCBI	NM_182919
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...