


Product Datasheet

ALIX (PDCD6IP) antibody (orb1333275)

Description	Mouse monoclonal antibody to ALIX (PDCD6IP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PDCD6IP(NP_037506) produced in HEK293T cell.
Target	PDCD6IP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	95.8 kDa
Uniprot ID	Q8WUM4
NCBI	NM_013374
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...