
Product Datasheet

Hydroxysteroid (17 beta) Dehydrogenase 4 (HSD17B4) antibody (orb1333390)

Description

Mouse monoclonal antibody to Hydroxysteroid


Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IF, IHC, WB

Western blot analysis of HEK293T cells w...

Immunogen

Full length human recombinant protein of human HSD17B4(NP_000405) produced in HEK293T cell.

Target

HSD17B4

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype

IgG1

Clonality

Monoclonal

MW

79.5 kDa

Uniprot ID
[P51659](#)
NCBI
[NM_000414](#)
Dilution Range

WB 1:4000, IHC 1:150, IF 1:100

Expiration Date

12 months from date of receipt.