


Product Datasheet


CD163 antibody (orb1333718)

Description	Mouse monoclonal antibody to CD163
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CD163(NP_004235) produced in HEK293T cell.
Target	CD163
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	125.3 kDa
Uniprot ID	Q86VB7
NCBI	NM_004244
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.




HEK293T cells were transfected with the

...



Western blot analysis of extracts (35ug)...



Immunohistochemical staining of paraffin...