


Product Datasheet

IL36 gamma (IL36G) antibody (orb1333929)

Description	Mouse monoclonal antibody to IL36 gamma (IL36G)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human IL1F9(NP_062564) produced in HEK293T cell.
Target	IL36G
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	18.5 kDa
Uniprot ID	Q9NZH8
NCBI	NM_019618
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...