

Product Datasheet

Solo (SESTD1) antibody (orb1333978)

Description	Mouse monoclonal antibody to Solo (SESTD1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human SESTD1(NP_835224) produced in HEK293T cell.
Target	SESTD1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	79.2 kDa
Uniprot ID	Q86VW0
NCBI	NM_178123
Dilution Range	WB 1:200~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.

