


Product Datasheet

MEK3 (MAP2K3) antibody (orb1334018)

Description	Mouse monoclonal antibody to MEK3 (MAP2K3)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MAP2K3(NP_659731) produced in HEK293T cell.
Target	MAP2K3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	39.1 kDa
Uniprot ID	P46734
NCBI	NM_145109
Dilution Range	WB 1:200~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...