


Product Datasheet

NANP antibody (orb1335031)

Description	Mouse monoclonal antibody to NANP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human NANP(NP_689880) produced in HEK293T cell.
Target	NANP
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	27.6 kDa
Uniprot ID	Q8TBE9
NCBI	NM_152667
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...