

Product Datasheet

Cytochrome P450 2B6 (CYP2B6) antibody (orb1335193)

Description Mouse monoclonal antibody to Cytochrome P450 2B6 (CYP2B6)

Species/Host Mouse

Reactivity Canine, Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human CYP2B6(NP_000758) produced in HEK293T cell.

Target CYP2B6

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 56.1 kDa

Uniprot ID [P20813](#)

NCBI [NM_000767](#)

Dilution Range WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...