


Product Datasheet

Coatomer subunit delta (ARCN1) antibody (orb1335203)

Description	Mouse monoclonal antibody to Coatomer subunit delta (ARCN1)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ARCN1(NP_001646) produced in HEK293T cell.
Target	ARCN1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	57 kDa
Uniprot ID	P48444
NCBI	NM_001655
Dilution Range	WB 1:500~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...