

Product Datasheet

Cytochrome P450 17A1 (CYP17A1) antibody (orb1335881)

Description Mouse monoclonal antibody to Cytochrome P450 17A1 (CYP17A1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human CYP17A1(NP_000093) produced in HEK293T cell.

Target CYP17A1

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 57.2 kDa

Uniprot ID [P05093](#)

NCBI [NM_000102](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...