


Product Datasheet

beta IV Tubulin (TUBB4A) antibody (orb1336094)

Description	Mouse monoclonal antibody to beta IV Tubulin (TUBB4A)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TUBB4 (NP_006078) produced in HEK293T cell.
Target	TUBB4A
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	49.4 kDa
Uniprot ID	P04350
NCBI	NM_006087
Dilution Range	WB 1:500-2000, IHC 1:150, IF 1:100, FC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...