

Product Datasheet

Spermine synthase (SMS) antibody (orb1336120)

Description Mouse monoclonal antibody to Spermine synthase (SMS)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human SMS (NP_004586) produced in HEK293T cell.

Target SMS

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 41.1 kDa

Uniprot ID [P52788](#)

NCBI [NM_004595](#)

Dilution Range WB 1:500~2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...